

# Continuous Passive Motion

## CPM FLEXIEE 4.1

RCF1741 [Shoulder CPM]

your new life  
New technologies for medical devices



- FLEXIEE 4.1 [RCF1741] is the stable and durable product that is smaller and lighter than other shoulder joint rehabilitation exercise device and is convenient for small hospitals or rental business.
- Manipulation is unnecessary from the location of an arm with symmetric arm cradle.
- Manual height adjuster is equipped, but gas spring is also installed that safe manipulation is available.
- With the entire weight of 12kgf, this device can be easily moved without additional equipment.
- With simple setup, rehabilitation exercise is available.
- In case of defect of the device, warranty service is applied up to 2 years.

**CARE TECH**  
[www.care-tech.co.kr](http://www.care-tech.co.kr)

### CPM FLEXIEE 4.1 Specification

|                            |                                       |
|----------------------------|---------------------------------------|
| Shoulder Flexion           | 180°                                  |
| Shoulder Extension         | 30°                                   |
| Timer                      | 0~60min                               |
| Speed                      | 185~285°/min                          |
|                            | Normal                                |
| Mode                       | Vibration                             |
|                            | Progress                              |
|                            | Limited Pause                         |
| Height [Gas spring/Stroke] | 250mm                                 |
| Weight                     | 12kgf                                 |
| Dimension                  | 470x620x965cm                         |
| Rated voltage              | 100~240VAC                            |
| Rated frequency            | 50/60Hz                               |
| PowerConsumption           | 75VA                                  |
| Electric safety class      | Class 1, Type B                       |
| Safety test                | IEC60601-1:2012,<br>IEC60601-1-6:2013 |




### Key-operation features

FLEXIEE series has the same key-operation process. Handy controller has the general function keys. Control panel on the main body has the special function keys.



▶ Handy controller of the FLEXIEE series



▶ Control panel of the FLEXIEE series main body

### Available both arms rehabilitation



▶ Operation of the shoulder direction of movement is inconvenient FLEXIEE 4.1 is available rehabilitation in desired direction without additional operations

### exercise with elbows bent



▶ Elbow angle is adjustable with simple resetting of arm link.

**CARE TECH**

[www.care-tech.co.kr](http://www.care-tech.co.kr)

B102 Tech dong, Sknatechpark #124, sagimakgol-ro,  
Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea  
462-721 CARETECH, Co.,Ltd.  
Tel : 82-31-767-3072 FAX : 82-31-767-3078